



CPAP

Continuous

Positive

Airway

Pressure

... bringing change to life[®]



What is CPAP?

Continuous Positive Airway Pressure (CPAP) is a non-invasive method to provide respiratory support to certain patients. CPAP has been shown to rapidly improve vital signs, gas exchange, the work of breathing, decrease the sense of dyspnea, and decrease the need for endotracheal intubation in patients who suffer from shortness of breath from Congestive Heart Failure (CHF) and Acute Pulmonary Edema (APE).

... bringing change to life®




Key Points of CPAP

- CPAP has been successfully demonstrated as an effective adjunct in the management of Pulmonary Edema secondary to CHF
- CPAP allows/buys time for administered medications to be able to work

... bringing change to life®

Lung Inflation


CPAP works by providing increased continuous gas pressures at the level of the *lower* airway structures, improving gas exchange in the alveoli



... bringing change to life®

CPAP Mechanism

- CPAP increases the airway pressures allowing for better gas diffusion & for re-expansion of collapsed alveoli
- CPAP allows the refilling of collapsed, airless alveoli
- CPAP expands the surface area of the collapsed alveoli allowing more surface area to be in contact with capillaries for gas exchange



... bringing change to life[®]



Goal of Therapy With CPAP

- Increase the amount of inspired oxygen
- Decrease the Work of Breathing (WOB)

In turn to:

- ❖ Decrease the need for intubation
- ❖ Decrease the hospital stay
- ❖ Decrease the mortality rate

... bringing change to life®



Contraindications

- O2-RESQ System
 - Pneumothorax
 - Decreased cardiac output and gastric distention
 - Severe facial injury (noninvasive use)
 - Respiratory Arrest
 - Hypotension secondary to Hypovolemia
- Always follow local established protocol

... bringing change to life®



The Requirements of CPAP

- The real requirement is for Continuous **CONSTANT** Positive Airway Pressure with **minimal** fluctuations in system pressure between inspiration and expiration
- The CONSTANT pressure will help reduce the patient WOB, which is important since it reduces the resistance to airflow that is present in lung disease

... bringing change to life®



O2-RESQ™ SYSTEM



... bringing change to life®



O2-RESQ SYSTEM

- Fixed Flow Generator with filter
- 72" Corrugated Anti-Asphyxia Circuit
- O2-CPAP™ Valve
- Mask & Head Strap
- Single Patient Use



... bringing change to life®



High Flow Generator

- O2-RESQ Generator is a fixed flow Venturi device that uses an oxygen supply in conjunction with entrained air to generate an output flow
 - Venturi tube ratio of approx. 10:1
 - Flow up to 140 Lpm and 30% FiO₂
- O2-RESQ Generator with DISS fitting hooks up to a 50psi Oxygen source
- Built in filter to protect device



... bringing change to life®

O₂-CPAP Valve & Anti-Asphyxia Valve

- O₂-CPAP Valves provide fixed pressure at any flow rate
 - Ensure that valve remains open during inspiration
- O₂-CPAP Valves are attached onto the elbow of the circuit
- 5.0, 7.5 & 10.0cm H₂O are the most commonly used pressures
- O₂-CPAP Valves can be easily changed to another pressure
 - 2.5, 5.0, 7.5, 10.0, 12.5, 15.0, 20.0cm H₂O
- Circuit has an Anti-Asphyxia valve which allows fresh air intake and prevents exhalation into the tubing in the event of no fresh gas flow



... bringing change to life[®]



O2-RESQ Face Mask Options

... bringing change to life®



SealMask™ and Head Strap

- SealMask
 - Comes with a Standard Non Adjustable Foam Forehead Pad
 - Has a single lip ribbed thermoplastic cushion for a comfortable low pressure seal
 - Three sizes available:
Large, Medium, & Small
- Head Strap
 - Attachment for Mask Posts
 - Dual Positioning Tabs



... bringing change to life®



AirMask™ and Head Strap

- AirMask
 - Comes with a Standard Non Adjustable Foam Forehead Pad
 - Utilizes a BVM style bladder cushion for a comfortable low pressure seal
 - Three sizes available:
Large, Medium, & Small
- Head Strap
 - Attachment for Mask Posts
 - Dual Positioning Tabs



... bringing change to life®



SoftMask™ and Head Strap

- SoftMask
 - Comes with a Standard Non Adjustable Foam Forehead Pad
 - Has a supple thermoplastic bladder cushion with valve for a comfortable low pressure adjustable seal
 - Three sizes available: Large, Medium, & Small
- Head Strap
 - Attachment for Mask Posts
 - Dual Positioning Tabs



... bringing change to life®



BiTrac ED™ Mask and Head Strap

- BiTrac ED™ Mask
 - Multi-positioning (up/down , in/out) OmniClip™ with Silicone Forehead Pad
 - Dual lip silicone cushion for excellent and comfortable low pressure seal
 - Four sizes available: X-Large, Large, Medium, & Small
- Head Strap
 - Attachment Clip to Mask
 - Dual Positioning Tabs



... bringing change to life®



BiTrac Shield™ ED Mask

- BiTrac Shield ED™ Mask
 - Single lip silicone face plate for comfortable low pressure seal
 - Excellent for patients with beards and facial issues
 - One size available: Adult
- Head Strap
 - Attachment Clip to Mask
 - Dual Positioning Tabs
 - Has a strap going the over the head to better position the mask



... bringing change to life®



O2-RESQ Accessories

... bringing change to life®

O2-RESQ CPAP Valves

- O2-CPAP Valves provide fixed pressure at any flow rate
- O2-CPAP Valves are attached onto the elbow of the circuit
- 5.0, 7.5 & 10.0cm H₂O are the most commonly used pressures
- O2-CPAP Valves can be easily changed to another pressure
 - Available pressures are:
 - 2.5, 5.0, 7.5, 10.0, 12.5, 15.0, 20.0cm H₂O



Always ensure that the valve remains open during inspiration

... bringing change to life®



Valved "T" Adapter

The Valved "T" is an optional accessory used to attach a nebulizer into the O2-RESQ



... bringing change to life®



O2-RESQ Filter


- The O2-RESQ Filter is an optional accessory that reduces the risk of medic infections by filtering exhaled air from patient
 - Bacterial/Viral
Efficiency 99.9999%
- Attach the Filter between the Circuit Elbow and the Mask



... bringing change to life®



O2-RESQ Filter Position



... bringing change to life[®]



O₂-RESQ Oxygen Hook Up

... bringing change to life[®]



O2-RESQ Connection to a Regulator



... bringing change to life®



Note: Regulator must be set to “0”



... bringing change to life®



Typical Flow Meter Connection

-Set Flow to Flush-



... bringing change to life®



“D” Cylinder O₂ Depletion Test Results

O₂-RESQ vs. Boussignac

- The elapsed time for the O₂-RESQ to run a “D” cylinder completely empty when connected to the DISS fitting on a regulator is 28 minutes.
- The elapsed time for the Boussignac to run a “D” cylinder completely empty when connected to the nipple fitting on a regulator set to 25 LPM, is 14 minutes.

Note: See discussions in the following slides for more information

... bringing change to life®



O2-RESQ “D” Cylinder Depletion Test Results

The elapsed time for the O2-RESQ to run a “D” cylinder completely empty when connected to the DISS fitting on a regulator is 28 minutes.

Test set up:

A BiTrac ED non-vented full face mask was securely fastened to a plastic test fixture and checked for leakage. The mask was connected to the O2-RESQ 72” corrugated anti-asphyxia circuit and thus to the fixed flow O2-RESQ Generator. A 10cmH2O O2-RESQ CPAP valve was attached to the circuit elbow. CPAP levels were continuously measured via a pressure tap located on a plastic test fixture. Testing was conducted with the fixed Generator connected directly to the 50 psi DISS fitting on the oxygen “D” cylinder regulator. The supply valve was fully opened and timing was initiated. The initial CPAP level generated was 10.1 cmH2O. **The CPAP level was stable at 10 cmH2O (+/- 0.2) throughout the test until the cylinder ran completely empty.**

... bringing change to life®



Boussignac “D” Cylinder Depletion Test Results

The elapsed time for the Boussignac to run a “D” cylinder completely empty when connected to the nipple fitting on a regulator set to 25 LPM, is 14 minutes - * Boussignac information.

Note: Pulmodyne internal testing found that the initial CPAP level generated by the Boussignac product was only 8.2 cmH₂O and decreased proportionately as the cylinder contents pressure decreased – approximately 0.3 – 0.4 cmH₂O / 250 psi decrease in cylinder pressure. Measured CPAP pressure was approx. 6.3 cmH₂O during the last 5 minutes of testing.

... bringing change to life[®]



*Boussignac “D” Cylinder Depletion Test Results

Q: How long will a cylinder of O₂ last using the Boussignac CPAP System?

A: Based on the CPAP level administered (8.5 to 10.0 cmH₂O@ 25 lpm) and the corresponding flow rate required using the Boussignac CPAP System, the following chart provides the approximate time that a full cylinder will last:

<u>Minutes of Oxygen by Cylinder Size</u>		All based on full 2200 PSI Cylinders	
Flow (LPM)	D Cylinder EMS Portable	E Cylinder EMS Portable	MCylinder EMS Ambulances
5	70	123	703
6	58	102	598
8	44	77	498
10	35	61	374
12	29	51	299
15	23	41	199
20	16	29	175
25	<u>14</u>	23	140

This information is taken directly from the Boussignac FAQ's found on www.vitaid.com

... bringing change to life®



Setting up the O2-RESQ System

Step 1 - Connect Generator directly to a 50psi Oxygen source on a Flow Meter or to a DISS fitting on an O2 tank regulator

Step 2



Remove the mask from Inner bag and attach the circuit to the mask. Turn on oxygen. Hand mask to patient for them to place on their face to get used to the flow. Then proceed to put the Head Strap on the patient.

Step 3



For all mask styles, Flip Head Strap forward and place mask on the patient's face. If using the BiTrac ED Mask, pinch the OmniClip, *slide it up or down* to find the best position on patient's forehead.

Step 4



For all mask styles, Flip Head Strap back over the patient's head, bring tabs forward on the top Head Strap and adjust equally to proper fit. Fold straps back to attach. Next attach the bottom 2 clips and repeat above sequence.

Step 5



Finally if using the BiTrac ED Mask, adjust the OmniClip *in and out* on Head Strap and Mask for best fit: Do Not Over Tighten the Head Strap. Monitor the Patient.

... bringing change to life®



Coaching The Patient

- Success is highly dependent upon patients tolerance and medics ability to coach the patient
 - Explain the procedure to the patient
 - Anticipate and control anxiety
 - CPAP may produce anxiety in some patients
 - Verbally coach breathing as needed
 - Consider having the patient hold the mask in place for a minute or so to reduce anxiety. As an option, the medic may hold the mask in place if a good seal is obtained.
 - Attach Head Straps loosely at first and gradually tighten until air leaks are eliminated
- Continue to coach patient

... bringing change to life[®]



CPAP vs. Intubation

- CPAP
 - Non-invasive
 - Easily discontinued
 - Easily adjusted
 - Use by EMT-B
 - Does not require sedation
 - Comfortable
- Intubation
 - Invasive
 - Usually don't extubate in field
 - Potential for infection
 - Requires highly trained personnel
 - Requires Sedation
 - Traumatic

... bringing change to life®